


## ScopeCheck - protein test

The Valisafe ScopeCheck protein test is a rapid test that has the ability to detect residual proteins left behind on the surfaces of washer-disinfectors, ultrasonic cleaners, endoscopes and other hard to clean surgical instruments. This protein test is based on a dye-binding solution used in clinical chemistry and can detect protein residues within 1 µg sensitivity.

- Conforms to BS EN ISO 15883
- Result within 10 seconds
- Easy to use - no incubation required
- Clear colour change
- Cost effective
- For use with:
  - Endoscopes
  - WD / Ultrasonic surfaces
  - Surgical instruments

<b>Name:</b>	Valisafe ScopeCheck
<b>Conformity:</b>	BS EN ISO 15883, HTM CFPP0106
<b>Time:</b>	Result within 10 Seconds
<b>Storage conditions:</b>	Refrigerated (2-8C)
<b>Lot no:</b>	Yes
<b>Correct symbols:</b>	Yes, in accordance with the standard
<b>Indicator:</b>	Diagnostic dye-binding vial
<b>Expiry date:</b>	6 months non-refrigerated 24 months refrigerated
<b>Packaging:</b>	25 tests per box including swabs plus extra, longer swabs if required. Plus 2 control tests



To determine if protein is present on a surface, swab the surface, place the swab in the colour changing solution and swirl the swab for 5 seconds. The swab can then be discarded. If protein is present, the reagent will turn blue within 10 seconds. The colour change provides a semi-quantitative measure of surface cleanliness; the higher the level of contamination of protein (from sources such as bioburden), the darker and faster the colour changes to blue.

Order details:

V3505005	V3504095 + 25 single 2.8mm swabs
V3505010	V3504095 + 25 single 1.9mm swabs
V3505013	V3504095 + 25 single 3.7mm swabs