



- **NOCOTEST®** are test sticks designed to semi-quantitative highlight of hydrogen peroxide (H₂O₂) and of peroxide in solutions. It especially enables to check that the diffusion with an automated surface disinfection device is homogeneous.
- **NOCOTEST®** sticks contain an organic redox indicator reacting with hydrogen peroxide to form a compound of blue oxidations.

REFERENCES AND PACKAGING

Reference	Packaging
5000.010	100 sticks tube

STORAGE

- Avoid exposing the sticks to sunlight and moisture.
- Store container in a cool and dry place (under 30°C).
- Do not leave opened boxes in a room during a treatment.
- Shelf-life : In the closed original packaging : 2 years from manufacturing date.
Once opened : 3 months from opening date.

INSTRUCTIONS FOR USE

- Open sticks tube NOCOTEST®.
- Take only needed quantity of NOCOTEST®.
- Replace the cap immediately.
- Do not touch the test zone.
- Briefly immerse (about 1 second) reading head of the stick in water.
- Set the sticks at desired locations in the room to be treated.
- At the end of the disinfection process (diffusion + contact time), collect the sticks.
- When hydrogen peroxide is present, test zone of the stick will color.
- Compare color of the test zone with the color scale on the label.
- The value corresponds to color on color scale that is the closest from the one on the reading head.
- After a short time, the reduction mechanism of reagent begins, turning the indicator to light/dark brown.

IMPORTANT :

- Elimination: sticks can be thrown with household waste.

OXY'PHARM

829 rue Marcel Paul
94500 Champigny-sur-Marne
commercial@oxypharm.net
T : +33 1 45 18 78 70

www.oxypharm.net

